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INTENDED USE

The dfu EASY TEST ASFV is an immunoassay kit for rapid diagnosis, qualitative detection of African swine fever virus (ASFV) antigen in whole blood of swine.

PRINCIPLE

dfu EASY TEST ASFV is a chromatographic immunoassay kit for the qualitative detection of ASFV antigen in whole blood of swine. You can drop the sample into the sample well and then add the assay solution. During test, if ASFV antigens are present in the swine whole blood, these antigens will bind to ASFV specific antibody-gold particle conjugates and move on the membrane by capillary forces, and then shows a red line on the test line (T) due to the binding with ASFV specific antibodies which are already coated on the membrane. Appearance of this red line indicates a positive result, while their absence indicates a negative result. A red line will always appear at the control line (C) to serve as a suitable performance control. The test device can quickly and simply detect the ASFV antigen within 15 minutes after dropping the sample.

KIT COMPONENTS

Reagent	20 tests/kit
ASFV Ag device	20 devices
Assay solution (0.4 mL/vial)	20 tubes
Dropper	20 ea
Disposable 20 μ L capillary tube	20 ea
Instruction for use	1 sheet

PRECAUTIONS


- For *in-vitro* veterinary diagnostic use only.
- Read this instruction for use thoroughly and follow all steps strictly for successful use of the product.
- Make sure to read the results on time.
- Do not use test device and reagents after expiration date, discard the assay solution if it is contaminated with bacteria or mold.
- The test device is sensitive to humidity and heat so it should be stored in a sealed pouch until ready for use. Use the test device immediately after opening pouch.
- Wear personal protective equipment such as lab coat, goggles, and disposable gloves while performing the test. Wash hands thoroughly afterwards.
- Handle all specimens, test device safely and in accordance with relevant regulations because of the risk of infection.
- After using, all wastes should be sterilized with high-pressure steam at 121°C for ≥ 15 minutes or comparable methods.
- Do not smoke, eat or drink in areas where specimens or kit components are handled.
- Never use reagents from another kit.

SPECIMEN COLLECTION AND STORAGE

1. Use whole blood of swine as sample.
2. Whole blood: Collect the blood into the tube contained anticoagulant such as EDTA, heparin, citrate according to clinical laboratory standard.
3. The anticoagulated whole blood samples should be tested within 24 hours after collection. If the sample is not tested immediately, it should be stored at 2 - 8°C within 72 hours, do not freeze.

TEST PROCEDURE

1. Place all components and samples at room temperature within 15 - 30 minutes before test.
2. Release the test device from the aluminum pouch and place it on a clean, dry, flat surface.
3. Using disposable capillary tube to add 20 μ L of whole blood (take the sample to the mark) into the sample well.
4. After 1 - 3 minutes, using dropper add 3 drops (approx. 100 μ L) of the assay solution into the sample well.
5. Read the result at 15 minutes. Do not read the result after 15 minutes.

**INTERPRETATION OF RESULTS**

1. **Positive:** Two red lines are appeared at the test line (T) and control line (C) in the result window.

2. **Negative:** One red line is appeared at the control line (C) in the result window.

3. **Invalid:** The red line does not appear at the control line (C) in the result window.

- The control line (C) is necessary in checking whether the reaction occurs abnormally or not. Therefore, if the control line does not appear, repeat the test with another test device.

STORAGE AND EXPIRATION

- Store the test kit at 2 - 30°C, do not freeze or heat.
- Expiration date: 24 months from manufacturing date.

LIMITATION OF THE TEST

The dfu EASY TEST ASFV is manufactured for diagnostic and screening detection of ASFV antigen in whole blood of swine, so do not completely exclude the possibility of false positive or false negative result. The final result needs to be based on clinical symptoms or higher assay techniques and the veterinarian's experience.