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Điện thoại: (+84) 28 71098734

Website : www.dfu.com.vn

dfu Amies Transport Medium w/o Charcoal (2S)

IVD kit for sampling and bacterial transport medium

INTENDED USE

To collect and transport mucous membranes from the oral cavity and nasal cavity, throat, vagina, wounds for bacteriological examination.

GENERAL INFORMATION

dfu Amies Transport Medium w/o Charcoal (2S) is a medium for the safe transport of bacteria in samples to the laboratory while maintaining their viability. The 0.3% sodium chloride, potassium, calcium and magnesium salts in the medium help maintain bacterial permeability. Phosphate buffer maintains the pH of the medium. Sodium thioglycolate and 0.9% agar maintain a stable environment for bacteria.

COMPONENT OF IVD KIT

Tube of 4 mL medium (containing 2 swabs) 100 set/box
Instruction for use 1 sheet/box

MEDIUM COMPOSITION

Component	g/L
Sodium chloride	3.0
Sodium hydrogen phosphate	1.15
Potassium dihydrogen phosphate	0.2
Potassium chloride	0.2
Sodium thioglycolate	1.0
Calcium chloride	0.1
Magnesium chloride	0.1
Agar	9.0
(pH 7.2±0.2)	

PRECAUTIONS

- For in vitro diagnostic use only.
- Observe aseptic techniques when using the product and transport promptly to the laboratory after completing sample collection.
- The specimen must be maintained at room temperature until testing.
- Do not break the tube.
- Do not use if the packaging is damaged or opened.
- Stop using if the reagent color has changed before use.
- The specimen may contain pathogenic microorganisms, so be careful when handling it.
- After use, sterilize the waste media and dispose of it as infectious agents.
- This product is disposable, do not reuse.

INSTRUCTION FOR USE

- Preparation
 - Specimens must be protected from contamination.
 - Prepare a sterilized specimen collection device and container.
- Test procedure
 - Tear off the plastic pouch ► direction until the tube cap is visible.
 - Remove the cap from transport medium tube.
 - Collect the specimen using the swab provided with the tube.
 - Immediately insert the swab with the collected specimen into the transport medium tube and close the cap.
 - Label the tube with patient name, specimen type, etc., and send it promptly to the laboratory.

3. Result evaluation

- After culturing the transported specimen, read the results.

QUALITY CONTROL

Control strain	Result
<i>Staphylococcus aureus</i> ATCC 25923	Grows well on culture medium
<i>Escherichia coli</i> ATCC 25922	Grows well on culture medium

1. Preparation for Bacterial Strains

- Prepare a standard strain at a concentration of 10^8 CFU/mL by comparing it with the turbidity of a 0.5 McFarland.
- Then dilute the standard strain in saline to a concentration of 10^6 CFU/mL.

2. Test procedure

- Dip a sterile swab into the diluted bacterial suspension, then place the swab in a transport medium tube and close the cap.
- After 24 hours, streak the swab onto Tryptone Soya Agar (TSA) and incubate at 35–37°C for 18–24 hours and examine the results.

STORAGE AND EXPIRATION

- dfu Amies Transport Medium w/o Charcoal (2S)** should be stored at 2 – 25°C.
- The expiration date is 1 year after its manufacturing date.

STERILIZATION INFORMATION

- The product is supplied in a sterile condition. Sterilization method: E-beam.
- Handling of Damaged Sterile Packaging:
 - If the sterile packaging is found to be torn, unsealed, damp, or shows any signs of compromised sterility before use, do not use the product.
 - Dispose of the product in accordance with medical waste management regulations.

SYMBOL EXPLANATION

 IVD	In vitro diagnostic medical device	 100	Sufficient for 100 sets
 REF	Catalogue / Product code		Consult instruction for use
 2	Do not reuse	 2°C – 25°C	Storage condition at 2 – 25°C

The owner and manufacturer:

DIAGNOSIS F CO., LTD

Address: 5th Floor, Factory No.02, Lot I-3b-1, N6 Street,
Saigon Hi-Tech Park, Tan Phu Ward, Thu Duc City, HCMC, Viet Nam

Product code: H52

MB170-Ver.01 (Rev.15/05/2025)

Tel: (+84) 28 71098734

Website : www.dfu.com.vn